

HabiStat

Leopard Gecko

Starter Kit



HSLGSKI

GB Leopard Gecko Starter Kit Assembly Instructions:

1: Insert the Heat Mat

- Slide the Heat Mat between the two panes of glass at the back of the Terrarium.
- Do not remove the foam insulation from in between the two glass panes. This will help insulate and ensure the heat is directed into the Terrarium. Please read the Heat Mat instructions carefully for any additional important information.
- Plug the Heat Mat into the socket of the Mat Stat (Thermostat).

2: Position the Mat Stat Probe

- Place the Mat Stat probe into the vent cable recess. You will need to remove a tab to fit the probe in.
- Ensure the probe is positioned over the Heat Mat and above the substrate. This will help regulate the surface temperature.
- The probe should be placed towards one side of the Terrarium to create a temperature gradient within the enclosure

3: Add Substrate

- Cover the bottom of the Terrarium with the substrate.
- Ensure there is only a thin layer of substrate over the position of the heat mat. This helps to avoid excessive heating and keeps the temperature in check.
- Please read the Heat Mat instructions carefully for any additional important information.

4: Attach the Mini UVB Kit

- Attach the Mini UVB kit to the inside lid of the tank using the provided clips.
- Place the UVB kit cable in one of the vent cable recesses after removing a tab.

5: Monitor Temperature

- Plug the Mat Stat into the mains power outlet and set the desired temperature using the adjustment knob.
- Use the included thermometer provided in the kit to monitor the temperature inside the Terrarium.
- Fix the thermometer next to the thermostat probe.
- Adjust the Mat Stat setting as necessary until the desired temperature is reached and maintained.

Note: The products enclosed in the kit may vary based on availability, but the basic setup and instructions should remain the same. Please follow all safety guidelines and read the individual instructions provided for each product for proper usage and care of your leopard gecko.


CONTENTS



Terrarium
61 x 38 x 30cm
(24 x 15 x 12")



Mat Stat 300 Watt



Heat Mat
28 x 28cm
(11 x 11") 12 Watt



Mini UVB Kit, 7%UVB
ShadeDweller Lamp, 8W



Sphagnum Moss



Leopard Gecko
Bedding 5 kg



Plastic Hide



Water Bowl



Bactericidal Cleaner



Coconut Cave x 2



Calcium Dish



Thermometer



Pure Calcium



Repton



Insecta Snack



Feeding Tongs